

Anti-Znf766 Rabbit Polyclonal Antibody For Wb And Elisa - Q5Hy98

Item Number: CM0024639

Introduction

ZNF766 rabbit polyclonal antibody validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Targets a 55kDa nuclear zinc finger protein implicated in transcriptional regulation. Ideal for gene expression research.

[Learn More](#)

Parameter	Value
Product Name	ZNF766 Rabbit pAb
Remarks/Alias	ZNF766
Species	Human
GeneID (Human)	90321
GeneID	90321
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-190 of human ZNF766 (NP_001010851.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q5HY98
Protein Weight	55kDa
Shipping	Ice bag