

Anti-Znf768 Polyclonal Antibody For Wb, Elisa - Q9H5H4

Item Number: CM0023356

Introduction

Rabbit polyclonal antibody targeting human ZNF768, suitable for Western blot and ELISA. Detects mouse ZNF768. Ideal for chromatin binding and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF768 Rabbit pAb
Remarks/Alias	ZNF768
Species	Human
GeneID (Human)	79724
GeneID	79724
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-270 of human ZNF768 (NP_078947.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9H5H4
Protein Weight	60kDa
Shipping	Ice bag