

Anti-Znf777 Polyclonal Antibody For Wb, Elisa - Q9Uld5

Item Number: CM0023207

Introduction

ZNF777 Rabbit pAb validated for WB and ELISA with cross-reactivity to human and mouse. Ideal for research on transcriptional repression and cell proliferation inhibition.

[Learn More](#)

Parameter	Value
Product Name	ZNF777 Rabbit pAb
Remarks/Alias	ZNF777
Species	Human
Gene ID	27153
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 532-831 of human ZNF777 (NP_056509.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9ULD5
Protein Weight	94kDa
Shipping	Ice bag