


Anti-Zo-1 Rabbit Monoclonal Antibody For If/Icc, Ihc-P, Elisa - Q07157

Item Number: CM0011398



Introduction

Rabbit monoclonal ZO-1 antibody for IF/ICC, IHC-P, ELISA. Recognizes human, mouse, rat ZO-1 (Q07157). Essential for tight junction research and epithelial barrier studies.

[Learn More](#)

Parameter	Value
Product Name	ZO-1 Rabbit mAb
Remarks/Alias	ZO-1
Species	Human
GeneID (Human)	7082
GeneID	7082
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human ZO-1 (NP_009093.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07157
Protein Weight	195kDa
Shipping	Ice bag