


Anti-Zo-1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q07157

Item Number: CM0017315



Introduction

High-quality rabbit polyclonal antibody against ZO-1 (Tight junction protein 1). Validated for Western blot, immunofluorescence-paraffin, and ELISA applications. Cross-reacts with human, mouse, rat. Ideal for tight junction research.

[Learn More](#)

Parameter	Specification
Product Name	ZO-1 Rabbit pAb
Remarks/Alias	ZO-1
Species	Human
Gene ID (Human)	7082
Gene ID	7082
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 18-273 of human ZO-1 (NP_003248.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07157
Protein Weight	195kDa
Shipping	Ice bag