

Anti-Zp2 Polyclonal Antibody For Wb - Q05996

Item Number: CM0014475



Introduction

Anti-ZP2 Polyclonal Antibody generated in rabbit, validated for WB, IF-P, and ELISA. Cross-reacts with human, mouse, rat ZP2. Use for oocyte zona pellucida and fertilization studies.

[Learn More](#)

Parameter	Value
Product Name	ZP2 Rabbit pAb
Remarks/Alias	ZPA; Zp-2; OOMD6; OZEMA6; ZP2
Species	Human
Gene ID	7783
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 651-745 of human ZP2 (NP_003451.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05996
Protein Weight	82kDa
Shipping	Ice bag