

Anti-Zbpw Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Bs86

Item Number: CM0022453

Introduction

High-quality rabbit polyclonal antibody against human ZBPW, validated for WB, IF-P, IHC-P, ELISA. Suitable for detecting ZBPW in human, mouse, rat samples. Involved in acrosome compaction and sperm-oocyte interaction.

[Learn More](#)

Parameter	Value
Product Name	ZBPW Rabbit pAb
Remarks/Alias	ZBPW1; SPGF66; ZBPW
Species	Human
Gene ID (Human)	11055
Gene ID	11055
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 45-200 of human ZBPW (NP_008940.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BS86
Protein Weight	40kDa
Shipping	Ice bag