

Anti-Zscan1 Polyclonal Antibody For Wb And Elisa - Q8Nbb4

Item Number: CM0021802

Introduction

Rabbit polyclonal antibody targeting human ZSCAN1 (also known as MZF-1, ZNF915). Validated for western blot and ELISA, recognizing a nuclear transcription factor of ~45 kDa. Suitable for studying transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	ZSCAN1 Rabbit pAb
Remarks/Alias	MZF-1; ZNF915; ZSCAN1
Species	Human
GeneID	284312
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 300-400 of human ZSCAN1 (NP_872378.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8NBB4
Protein Weight	45kDa
Shipping	Ice bag