

Anti-Zscan21 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y5A6

Item Number: CM0019076

Introduction

Rabbit polyclonal antibody targeting human ZSCAN21 (ZNF38, Zipro1), validated in WB and ELISA. Cross-reacts with mouse and rat. Suitable for gene regulation and spermatogenesis research.

[Learn More](#)

Parameter	Value
Product Name	ZSCAN21 Rabbit pAb
Remarks/Alias	ZNF38; Zipro1; NY-REN-21; ZSCAN21
Species	Human
Gene ID (Human)	7589
Gene ID	7589
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-280 of human ZSCAN21 (NP_666019.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5A6
Protein Weight	54kDa
Shipping	Ice bag