

Anti-Zswim1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Br11

Item Number: CM0013445

Introduction

Rabbit polyclonal antibody targeting human ZSWIM1, validated for Western blot and ELISA. Recognizes mouse ZSWIM1. For studying zinc finger SWIM domain proteins and diagnostic reagent development.

[Learn More](#)

Parameter	Value
Product Name	ZSWIM1 Rabbit pAb
Remarks/Alias	C20orf162; ZSWIM1
Species	Human
Gene ID (Human)	90204
Gene ID	90204
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 411-485 of human ZSWIM1 (NP_542170.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9BR11
Protein Weight	55kDa
Shipping	Ice bag