

Anti-Zwint Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O95229

Item Number: CM0030543



Introduction

Rabbit polyclonal antibody targeting human ZWINT, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, rat. Detects outer kinetochore KNL1 complex subunit ZWINT, essential for chromosome segregation. Ideal for mitosis and cell cycle studies.

[Learn More](#)

Parameter	Specification
Product Name	ZWINT Rabbit pAb
Remarks/Alias	SIP30; ZWINT1; KNTC2AP; HZwint-1; ZWINT
Species	Human
GeneID (Human)	11130
GeneID	11130
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-277 of human ZWINT (NP_008988.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95229
Protein Weight	31kDa
Shipping	Ice bag