

Anti-Zyxin Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15942

Item Number: CM0022358

Introduction

Rabbit polyclonal antibody to human Zyxin (ESP-2, HED-2), suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for focal adhesion and cytoskeleton studies.

[Learn More](#)

Parameter	Value
Product Name	Zyxin Rabbit pAb
Remarks/Alias	ESP-2, HED-2, Zyxin
Species	Human
Gene ID (Human)	7791
Gene ID	7791
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-572 of human Zyxin (NP_003452.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15942
Protein Weight	61kDa
Shipping	Ice bag