


Anti-Zzz3 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Iyh5

Item Number: CM0012894



Introduction

Rabbit polyclonal antibody against human ZZZ3 (ATAC1) for WB and ELISA. Cross-reacts with human and rat. Useful for studying histone acetylation and ATAC complex function in transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	ZZZ3 Rabbit pAb
Remarks/Alias	ATAC1; ZZZ3
Species	Human
Gene ID (Human)	26009
Gene ID	26009
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ZZZ3 (Q8IYH5)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8IYH5
Protein Weight	102kDa
Shipping	Ice bag