


Antxr2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P58335

Item Number: CM0025863



Introduction

ANTXR2 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat. Targets anthrax toxin receptor 2 (CMG2), suitable for studying anthrax pathogenesis and extracellular matrix interactions.

[Learn More](#)

Parameter	Specification
Product Name	ANTXR2 Rabbit pAb
Remarks/Alias	HFS; ISH; JHF; CMG2; CMG-2; ANTXR2
Species	Human
Gene ID (Human)	118429
Gene ID	118429
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-300 of human ANTXR2 (NP_001139266.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P58335
Protein Weight	54kDa
Shipping	Ice bag