


Ap2B1 Rabbit Monoclonal Antibody For Wb And Elisa - P63010

Item Number: CM0007354



Introduction

AP2B1 rabbit monoclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat. Targets the beta subunit of the AP-2 clathrin adaptor complex involved in endocytosis. Ideal for vesicle trafficking and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	AP2B1 Rabbit mAb
Remarks/Alias	AP2B1, ADTB2, AP105B, AP2-BETA, CLAPB1, AP-2 complex subunit beta
Species	Human
Gene ID	163
Immunogen	Recombinant fusion protein containing amino acids 752-951 of human AP2B1 (NP_001025177.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63010
Protein Weight	98kDa/104kDa/105kDa
Shipping	Ice bag