


Ap2B1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P63010

Item Number: CM0022064



Introduction

AP2B1 rabbit polyclonal antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Targets AP-2 complex subunit beta involved in clathrin-mediated endocytosis. Ideal for neuroscience and membrane trafficking research.

[Learn More](#)

Parameter	Value
Product Name	AP2B1 Rabbit pAb
Remarks/Alias	ADTB2; AP105B; CLAPB1; AP2-BETA; AP2B1
Species	Human
GeneID (Human)	163
GeneID	163
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 752-951 of human AP2B1 (NP_001025177.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63010
Protein Weight	105kDa
Shipping	Ice bag