

Ap3B1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O00203

Item Number: CM0027171



Introduction

AP3B1 rabbit polyclonal antibody validated for Western blot (WB), IHC-P, and ELISA.

Recognizes human, mouse, and rat AP-3 complex subunit beta-1 (O00203). Ideal for studying lysosomal protein sorting and related diseases.

[Learn More](#)

Parameter	Value
Product Name	AP3B1 Rabbit pAb
Remarks/Alias	PE; HPS; HPS2; ADTB3; ADTB3A; AP3B1
Species	Human
GeneID (human)	8546
GeneID	8546
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 895-1094 of human AP3B1 (NP_003655.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00203
Protein Weight	121kDa
Shipping	Ice bag