


Apaf1 Rabbit Pab For Wb, If/Icc, Elisa - O14727

Item Number: CM0018635



Introduction

APAF1 Rabbit pAb (O14727) is a polyclonal antibody for WB, IF/ICC, ELISA. Detects human, mouse, rat APAF1, a cytoplasmic protein mediating cytochrome c-dependent caspase-9 activation and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	APAF1 Rabbit pAb
Remarks/Alias	CED4; APAF-1; APAF1
Species	Human
GeneID (human)	317
GeneID	317
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 470-600 of human APAF1 (NP_863651.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14727
Protein Weight	38-142 kDa
Shipping	Ice bag