

Apc Anti-Cd141/Thrombomodulin Monoclonal Antibody For Flow Cytometry - P07204

Item Number: CM0005395

Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD141/Thrombomodulin, validated for flow cytometry. Immunogen covers extracellular domain (aa 19-515). Suitable for endothelial cell and coagulation research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD141/Thrombomodulin mAb
Remarks/Alias	TM; THRM; AHUS6; BDCA3; CD141; BDCA-3; THPH12
Species	Human
GeneID (human)	7056
GeneID	7056
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-515 of human CD141/Thrombomodulin (NP_000352.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07204
Protein Weight	60kDa
Shipping	Ice bag