


Apc Anti-Cd66B Monoclonal Antibody For Flow Cytometry - P31997

Item Number: CM0005866



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD66b (CEACAM8, NCA-95). The antibody recognizes the CD66b antigen, a cell adhesion molecule on neutrophils. It is provided as an APC conjugate for direct detection in flow cytometry applications.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD66b mAb
Remarks/Alias	CD67; CGM6; CD66b; NCA-95
Species	Human
GeneID (Human)	1088
GeneID	1088
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-319 of human CD66b (NP_001807.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P31997
Protein Weight	38kDa
Shipping	Ice bag