

Apc Anti-F4/80 Monoclonal Antibody For Flow Cytometry - Q61549

Item Number: CM0009126

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse F4/80 antigen. Validated for flow cytometry applications, suitable for detecting macrophages in murine samples. Swiss-Prot Q61549.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse F4/80 mAb
Remarks/Alias	Emr1; Ly71; F4/80; Gpf480; TM7LN3; DD7A5-7; EGF-TM7
Species	Mouse
GenelD (Human)	2015
GenelD	13733
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-644 of mouse F4/80 (NP_034260.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61549
Protein Weight	102 kDa
Shipping	Ice bag