

Apc Anti-Human Cd101/Igsf2 Monoclonal Antibody For Flow Cytometry - Q93033

Item Number: CM0011465

Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD101 (IGSF2, V7, EWI-101). Designed for flow cytometry, it detects human CD101, a T-cell proliferation inhibitor expressed on immune cells and dendritic cells.

[Learn More](#)

Parameter	Specification
Product Name	APC Rabbit anti-Human CD101/IGSF2 mAb
Remarks/Alias	V7; IGSF2; EWI-101
Species	Human
GenelD (Human)	9398
GenelD	9398
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-410 of human CD101/IGSF2 (NP_004249.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q93033
Protein Weight	115kDa
Shipping	Ice bag