


Apc Anti-Human Cd156C/Adam10 Monoclonal Antibody For Flow Cytometry - O14672

Item Number: CM0005537



Introduction

APC-conjugated rabbit monoclonal antibody targeting human ADAM10 (CD156c). Specifically designed for flow cytometry detection of ADAM10, a key sheddase involved in ectodomain shedding, Notch signaling, and Alzheimer's disease. Research-grade IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD156c/ADAM10 mAb
Remarks/Alias	RAK; kuz; AD10; AD18; MADM; CD156c; CDw156; HsT18717
Species	Human
Gene ID	102
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 214-672 of human ADAM10 (NP_001101.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O14672
Protein Weight	84kDa
Shipping	Ice bag