


Apc Anti-Human Cd163 Monoclonal Antibody For Flow Cytometry - Q86Vb7

Item Number: CM0010996



Introduction

APC-conjugated mouse monoclonal antibody targeting human CD163 (M130) scavenger receptor. Validated for flow cytometry detection of monocytes and macrophages. Ideal for immunophenotyping, inflammation, and macrophage activation studies.

[Learn More](#)

Parameter	Value
Product Name	APC Mouse anti-Human CD163 mAb
Remarks/Alias	M130; MM130; SCAR11
Species	Human
Gene ID	9332
Immunogen	Recombinant protein of human CD163
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q86VB7
Protein Weight	125kDa
Shipping	Ice bag