

Apc Anti-Human Cd24 Monoclonal Antibody For Flow Cytometry - P25063

Item Number: CM0005800

Introduction

APC-conjugated rabbit anti-human CD24 monoclonal antibody for flow cytometry. Targets the GPI-anchored glycoprotein CD24 (8-13 kDa), expressed on B-cells, involved in B-cell activation and differentiation.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD24 mAb
Remarks/Alias	CD24A; CD24
Species	Human
GeneID (Human)	100133941
GeneID	100133941
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-80 of human CD24 (NP_037362.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P25063
Protein Weight	8kDa; 13kDa
Shipping	Ice bag