


Apc Anti-Human Cd28 Rabbit Monoclonal Antibody For Flow Cytometry - P10747

Item Number: CM0008523



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD28 (Tp44), validated for flow cytometry. Ideal for studying T-cell activation, co-stimulation, and immune homeostasis. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD28 mAb
Remarks/Alias	Tp44
Species	Human
Gene ID (Human)	940
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P10747
Protein Weight	25kDa
Shipping	Ice bag