

Apc Anti-Human Cd335/Nkp46 Monoclonal Antibody For Flow Cytometry - O76036

Item Number: CM0005913



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD335/NKp46 (SWISS: O76036) specifically for flow cytometry. Reacts with human NK cells, useful for tumor cell lysis studies.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD335/NKp46 mAb
Remarks/Alias	LY94; CD335; NKP46; NK-p46
Species	Human
Gene ID (Human)	9437
Gene ID	9437
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-254 of human CD335/NKp46 (NP_004820.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O76036
Protein Weight	21kDa/22kDa/23kDa/32kDa/34kDa
Shipping	Ice bag