

Apc Anti-Human Cd79A Monoclonal Antibody For Flow Cytometry (Intracellular) - P11912

Item Number: CM0005583



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD79a (MB1, Ig-alpha). Designed for intracellular flow cytometry (FC) analysis of human B-cells. Recognizes the B-cell antigen receptor alpha chain. MW: 20/25 kDa.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD79a mAb
Remarks/Alias	IGA; MB1; MB-1; IGAAlpha
Species	Human
GeneID (human)	973
GeneID	973
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 165-225 of human CD79a (NP_001774.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P11912
Protein Weight	20kDa/25kDa
Shipping	Ice bag