


Apc Anti-Human Tlr3 Monoclonal Antibody For Flow Cytometry - O15455

Item Number: CM0011928



Introduction

APC-conjugated rabbit monoclonal antibody targeting human TLR3 (CD283, IIAE2, IMD83). Validated for intracellular flow cytometry, ideal for studying innate immune responses.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human TLR3 mAb
Remarks/Alias	CD283; IIAE2; IMD83
Species	Human
GeneID (Human)	7098
GeneID	7098
Immunogen	Recombinant protein corresponding to a sequence within amino acids 24-339 of human TLR3 (NP_003256.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	O15455
Protein Weight	104kDa
Shipping	Ice bag