

ApC Anti-Integrin Alpha 5 (Itga5/Cd49E) Monoclonal Antibody For Flow Cytometry - P08648

Item Number: CM0009991



Introduction

APC-conjugated rabbit monoclonal antibody targeting human/mouse Integrin alpha 5 (ITGA5/CD49e). Validated for flow cytometry. Recognizes fibronectin receptor involved in cell adhesion. Recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human/Mouse Integrin alpha 5 (ITGA5/CD49e) mAb
Remarks/Alias	FNRA; CD49e; VLA-5; VLA5A
Species	Human
GeneID	3678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-995 of human Integrin alpha 5 (ITGA5/CD49e) (NP_002196.4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	P08648
Protein Weight	115kDa
Shipping	Ice bag