


Apc Anti-Lgr5 Monoclonal Antibody For Flow Cytometry - 075473

Item Number: CM0010757



Introduction

APC-conjugated rabbit monoclonal antibody against human LGR5 (GPR49), a Wnt receptor and stem cell marker. Suitable for flow cytometry applications in intestinal epithelium, hair follicle, and oncology research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human LGR5 (GPR49) mAb
Remarks/Alias	FEX; HG38; GPR49; GPR67; GRP49
Species	Human
Gene ID (Human)	8549
Gene ID	8549
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 480-580 of human LGR5 (NP_003658.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	075473
Protein Weight	100kDa
Shipping	Ice bag