

Apc Anti-Mouse Cd11C Monoclonal Antibody For Flow Cytometry - Q9Qxh4

Item Number: CM0008554

Introduction

APC-conjugated anti-mouse CD11c monoclonal antibody from Armenian hamster. Validated for flow cytometry. Suitable for investigating monocyte adhesion, chemotaxis, and inflammatory response mediated by Integrin alpha-X.

[Learn More](#)

Parameter	Value
Product Name	APC Anti-Mouse CD11c mAb
Remarks/Alias	Cr4; N418; Cd11c
Species	Mouse
GeneID (Human)	3687
GeneID	16411
Immunogen	Recombinant protein
Source	Armenian Hamster
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QXH4
Protein Weight	129kDa
Shipping	Ice bag