

Apc Anti-Mouse Cd163 Monoclonal Antibody For Flow Cytometry - Q2Vlh6

Item Number: CM0012663

Introduction

This APC-conjugated rabbit monoclonal antibody against mouse CD163 (M130) is validated for flow cytometry. It recognizes the extracellular domain of mouse CD163, with a predicted molecular weight of 121-125 kDa.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD163 mAb
Remarks/Alias	CD163v2; CD163v3; M130
Species	Mouse
GenelD (human)	9332
GenelD	93671
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 39-260 of mouse CD163 (NP_444324.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q2VLH6
Protein Weight	121kDa/125kDa
Shipping	Ice bag