

Apc Anti-Mouse Cd19 Mab - Mouse Cd19

Item Number: CM0011116

Introduction

APC-conjugated rat monoclonal antibody targeting Mouse CD19 for flow cytometry. Specifically reacts with mouse CD19, enabling immune cell analysis in research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	APC Anti-Mouse CD19 mAb
Species	Mouse
Gene ID	12478
Immunogen	Recombinant protein of Mouse CD19
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
Shipping	Ice bag