

Apc Anti-Mouse Cd19 Monoclonal Antibody For Flow Cytometry - P25918

Item Number: CM0004495

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse CD19, validated for flow cytometry. Ideal for B-cell identification, differentiation, and signaling research. Recognizes extracellular domain, 60kDa.

[Learn More](#)

Parameter	Details
Product Name	APC Rabbit anti-Mouse CD19 mAb
Remarks/Alias	Cd19
Species	Mouse
GenelD (Human)	930
GenelD	12478
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-287 of mouse CD19 (NP_033974.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25918
Protein Weight	60kDa
Shipping	Ice bag