


ApC Anti-Mouse Cd44 Monoclonal Antibody For Fc - P15379

Item Number: CM0008878



Introduction

APC-conjugated rabbit anti-mouse CD44 monoclonal antibody validated for flow cytometry. Targets the hyaluronate receptor involved in lymphocyte homing, cell adhesion, and migration.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD44 mAb
Remarks/Alias	Ly-24; Pgp-1; HERMES
Species	Mouse
GeneID (Human)	960
GeneID	12505
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-176 of mouse CD44 (NP_033981.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P15379
Protein Weight	85kDa
Shipping	Ice bag