

Apc Anti-Mouse Cd49F Monoclonal Antibody For Fc - Q61739

Item Number: CM0008259

Introduction

Rat monoclonal antibody against mouse CD49f, conjugated to APC, suitable for flow cytometry (FC). Recognizes mouse and human integrin alpha-6 (Q61739). Ideal for studying cell adhesion and laminin receptor function.

[Learn More](#)

Parameter	Value
Product Name	APC Rat anti-Mouse CD49f mAb
Remarks/Alias	Cd49f; VLA-6; 5033401005Rik
Species	Mouse
GeneID (Human)	3655
GeneID	16403
Immunogen	Synthetic peptide derived from mouse CD49f
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Human
SWISS	Q61739
Protein Weight	122 kDa
Shipping	Ice bag