


APC-Conjugated Anti-Cd223/Lag-3 Rabbit Monoclonal Antibody For Flow Cytometry - P18627

Item Number: CM0012162



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD223/LAG-3 for flow cytometry. LAG-3 is an immune checkpoint receptor regulating T-cell activity. Suitable for immuno-oncology research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD223/LAG-3 mAb
Remarks/Alias	CD223
Species	Human
Gene ID (Human)	3902
Gene ID	3902
Immunogen	Recombinant protein: recombinant fusion protein containing a sequence corresponding to amino acids 23-430 of human LAG3 (NP_002277.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P18627
Protein Weight	57 kDa
Shipping	Ice bag