


Apc-Conjugated Anti-Siglec-7/Cd328 Monoclonal Antibody (Rabbit) For Flow Cytometry - Q9Y286

Item Number: CM0012594



Introduction

APC-conjugated rabbit monoclonal antibody targeting human Siglec-7 (CD328, AIRM1). Validated for flow cytometry (FC). Suitable for natural killer cell inhibition and immune receptor research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human Siglec-7/CD328 mAb
Remarks/Alias	p75; QA79; AIRM1; CD328; AIRM-1; CDw328; D-siglec; SIGLEC-7; SIGLECP2; SIGLEC19P; p75/AIRM1
Species	Human
GeneID (Human)	27036
GeneID	27036
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-144 of human Siglec-7/CD328 (NP_055200.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9Y286
Protein Weight	51kDa
Shipping	Ice bag