


Apc-Conjugated Rabbit Anti-Human Cd30 Monoclonal Antibody For Flow Cytometry - P28908

Item Number: CM0011903



Introduction

APC-conjugated rabbit monoclonal antibody against human CD30 (TNFRSF8, Ki-1 antigen). Suitable for flow cytometry analysis of CD30 expression in lymphoma and activated lymphocyte research.

[Learn More](#)

Product Name	APC Rabbit anti-Human CD30 mAb
Remarks/Alias	CD30; Ki-1; D1S166E
Species	Human
GeneID (Human)	943
GeneID	943
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 19-379 of human CD30 (NP_001234.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28908
Protein Weight	14kDa/50kDa/63kDa
Shipping	Ice bag