


Apc-Conjugated Rabbit Anti-Human Cd301/Clec10A Monoclonal Antibody For Flow Cytometry - Q8Iun9

Item Number: CM0012453



Introduction

APC-conjugated rabbit monoclonal antibody specific for human CD301/CLEC10A, validated for flow cytometry. Recognizes GalNAc-terminated glycans on dendritic cells, suitable for immunology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD301/CLEC10A mAb
Remarks/Alias	HML; MGL; HML2; CD301; CLECSF13; CLECSF14
Species	Human
GeneID (human)	10462
GeneID	10462
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 61-316 of human CD301/CLEC10A (NP_878910.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8IUN9
Protein Weight	28kDa/32kDa/35kD
Shipping	Ice bag