

APC-Conjugated Rabbit Anti-Human Cd3E Monoclonal Antibody For Intracellular Flow Cytometry - P07766

Item Number: CM0004151



Introduction

APC-conjugated rabbit monoclonal antibody targeting human CD3 epsilon chain (CD3E), validated for intracellular flow cytometry. Ideal for T-cell activation and immune response research. Reacts with human. Corresponds to UniProt P07766.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human cCD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
GeneID (Human)	916
GeneID	916
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag