

Apc-Conjugated Rabbit Anti-Human Cd62L/L-Selectin Monoclonal Antibody For Flow Cytometry - P14151

Item Number: CM0012805



Introduction

APC-conjugated rabbit monoclonal antibody against human CD62L/L-selectin (SELL). Designed for flow cytometry (FC) to detect L-selectin, a cell adhesion molecule mediating leukocyte tethering and rolling. Suitable for human samples.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD62L/L-Selectin mAb
Remarks/Alias	TQ1; LAM1; LEU8; LNHR; LSEL; CD62L; LYAM1; PLNHR; LECAM1
Species	Human
GeneID (Human)	6402
GeneID	6402
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-332 of human CD62L/L-Selectin (NP_000646.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P14151
Protein Weight	42kDa
Shipping	Ice bag