

Apc-Conjugated Rabbit Anti-Mouse Foxp3 Monoclonal Antibody For Flow Cytometry (Intracellular) - Q99Jb6

Item Number: CM0012414



Introduction

APC-conjugated rabbit monoclonal antibody against mouse Foxp3 (scurfin), validated for intracellular flow cytometry (FC intra). Detects mouse Foxp3 protein at ~47kDa. Ideal for studying regulatory T-cell (Treg) development and immune homeostasis.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse Foxp3 mAb
Remarks/Alias	sf; JM2; scurfen
Species	Mouse
GeneID (Human)	50943
GeneID (Mouse)	20371
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse Foxp3 (NP_473380.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	Q99JB6
Protein Weight	47kDa
Shipping	Ice bag