

## ApC Mouse Anti-Human Cd11B Monoclonal Antibody For Fc - P11215

Item Number: CM0011721



### Introduction

APC-conjugated mouse monoclonal antibody targeting human CD11b (ITGAM, MAC-1). Validated for flow cytometry, specific to human. Molecular weight 127 kDa. For immunology and innate immunity research.

[Learn More](#)

| Parameter        | Value                                  |
|------------------|----------------------------------------|
| Product Name     | APC Mouse anti-Human CD11b mAb         |
| Remarks/Alias    | CR3A; MO1A; CD11B; MAC-1; MAC1A; SLEB6 |
| Species          | Human                                  |
| Gene ID          | 3684                                   |
| Immunogen        | Recombinant protein of human CD11b     |
| Source           | Mouse                                  |
| Category         | Monoclonal Antibodies                  |
| Application      | FC                                     |
| Cross Reactivity | Human                                  |
| SWISS            | P11215                                 |
| Protein Weight   | 127kDa                                 |
| Shipping         | Ice bag                                |