

Apc Mouse Anti-Human Cd11B Monoclonal Antibody For Flow Cytometry - P11215

Item Number: CM0008220



Introduction

APC-conjugated mouse monoclonal antibody targeting human CD11b (Integrin alpha-M) validated for flow cytometry. Ideal for innate immunity, cell adhesion, and neutrophil/monocyte studies.

[Learn More](#)

Parameter	Value
Product Name	APC Mouse anti-Human CD11b mAb
Remarks/Alias	CR3A; MO1A; CD11B; MAC-1; MAC1A; SLEB6
Species	Human
GeneID (Human)	3684
GeneID	3684
Immunogen	Synthetic peptide / A synthesized peptide derived from human CD11b.
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P11215
Protein Weight	127kDa
Shipping	Ice bag