

ApC Mouse Anti-Human Cd11C Monoclonal Antibody For Flow Cytometry - P20702

Item Number: CM0010354



Introduction

APC-conjugated mouse monoclonal antibody targeting human CD11c (ITGAX) for flow cytometry. Reacts with human monocyte and granulocyte integrin alpha-X, suitable for inflammatory response and adhesion research.

[Learn More](#)

Parameter	Value
Product Name	APC Mouse anti-Human CD11c mAb
Remarks/Alias	CD11C; SLEB6
Species	Human
GeneID (Human)	3687
GeneID	3687
Immunogen	Recombinant protein of human CD11c.
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20702
Protein Weight	128kDa
Shipping	Ice bag