


Apc Mouse Anti-Human Cd1A Monoclonal Antibody For Flow Cytometry - P06126

Item Number: CM0012327



Introduction

APC Mouse anti-Human CD1a monoclonal antibody designed for flow cytometry. Detects CD1a, an antigen-presenting protein on T-cells, dendritic cells, and Langerhans cells. Reacts with Human. Supported by CamelBio comprehensive IVD raw material supply.

[Learn More](#)

Parameter	Value
Product Name	APC Mouse anti-Human CD1a mAb
Remarks/Alias	R4; T6; CD1; FCB6; HTA1
Species	Human
Gene ID (Human)	909
Gene ID	909
Immunogen	Synthetic peptide A synthesized peptide derived from human CD1a.
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06126
Protein Weight	37kDa
Shipping	Ice bag