

Apc Mouse Anti-Human Cd3 Monoclonal Antibody For Flow Cytometry - P07766

Item Number: CM0012148



Introduction

APC-conjugated mouse monoclonal antibody against human CD3 epsilon chain, validated for flow cytometry. Specific for T-cell surface receptor, with a predicted molecular weight of 23 kDa. Ideal for immune profiling, T-cell activation, and signaling studies.

[Learn More](#)

Specification	Details
Product Name	APC Mouse anti-Human CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
GeneID (Human)	916
GeneID	916
Immunogen	Recombinant protein Recombinant protein of human CD3.
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag