

APC Mouse Anti-Human CD42A (Gpix) Monoclonal Antibody For Flow Cytometry - P14770

Item Number: CM0003995



Introduction

APC-conjugated mouse monoclonal antibody targeting human CD42a (GPIX/GP9), validated for flow cytometry. Specific for the platelet glycoprotein IX subunit of the GPIb-V-IX complex. Suitable for hemostasis, platelet adhesion, and Bernard-Soulier syndrome research.

[Learn More](#)

Parameter	Specification
Product Name	APC Mouse anti-Human CD42a mAb
Remarks/Alias	GPIX; CD42a
Species	Human
GeneID (Human)	2815
GeneID	2815
Immunogen	Recombinant protein of human CD42a
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P14770
Protein Weight	19kDa
Shipping	Ice bag